



# JHEAR OTC hearing aid Platform introduction

---

**April. 2024**

**<http://en.jhearing.com>**


**[tech@jhearing.com](mailto:tech@jhearing.com)**

# Agenda

- JHAR one page introduction
- J10 platform introduction
- J10 Platform overview
- J10 OTC Hearing Aid Scenarios
- J10 Platform Performance picture (OTC)
- J10 design support resource
- J10 Hearing Aid Accessories
- J10 platform eco-system build for ramping up market
- J10 hearing aid platform advantages
- JHEAR Hearing aid roadmap
- J10 platform Success stories

# Wireless audio streaming

- Battery autonomy and audio BW landscape



# JHEAR corporate evolution



JHEAR

Focus

In 2004, established R&D team for VOIP and soft switch systems

In 2020, the intelligent Wireless Division was established. Provide a complete wireless communication solution, software system, terminal software, and IoT SaaS

In 2022, it launched the first hearing aid SiP chip J10 in China, integrating core functional chips such as DSP, BLE, MCU, etc., and won more than 10 patents and Software copyright




In 2015, the hearing aid department was established, to support onsemi platform business in APAC region


In 2024, 2<sup>nd</sup> generation J102 SiP tape out with more diversity applications

In 2010, data collection, mining and processing business. natural language text processing (NLP machine learning)

In 2021, the healthcare department was established, including hearing aid business and blood glucose business. The hearing aid business includes hardware, software, algorithms, fitting solution, also support CGM total solution as IDH.

# JHEAR business structure





Hardware

software

algorithm

Fitting tool


System design

# JHEAR J10 Smallest BLUETOOTH® 5.2 Wireless SiP




**JHEARING ELECTORICS** J10 BLUETOOTH® 5.2 Wireless Module is ultra-small, designed for connected portable devices in industrial and medical applications. This module offers ultra-low power consumption, integrated dual core- MCU core and DSP core for user applications. The J10 wireless module provides a **4mm x 3.5mm x 0.75mm LGA package**, 48MHz Arm® Cortex® M3, 50ohm match to external antenna, and various I/O interfaces such as UART, QSPI, SPI, GPIO, ADC, DAC, PWM, and I2C. The Wireless Module has an operating temperature range -40°C to +85°C and an input supply range of **1.2V to 3.6V**.


**J102 SiP module offer 4mm\*4mm\*0.75mm LGA package with more GPIOs**



## Hearing Aid Applications




<https://en.jhearing.com/>



# JHEAR J10 SiP & J102 SiP Package

## J10 SiP



- 4mm\*3.5mm\*0.75mm
- 20 pins
- Mass Production



Pin#	Pin Def	Pin Name	Pin Function	Pin Description
1	A1	TMS	I/O	Test Mode Select
2	A2	IO13	I/O	I/O
3	A3	IO14	I/O	I/O
4	A4	GND	GND	GND
5	A5	ANT	RF	RF signal input/output (Antenna)
6	B1	IO10	I/O	I/O
7	B2	TCK	I/O	External clock input
8	B3	IO11	I/O	I/O
9	B4	AGND	AGND	AGND
10	B5	AGND	AGND	AGND
11	C1	IO5	I/O	I/O
12	C2	IO3	I/O	I/O
13	C3	NRESET	I/O	Reset
14	C4	VDC	VDC	DC-DC output voltage to external LC filter
15	C5	VBAT	Power Source	Power Source
16	D1	IO1	I/O	I/O
17	D2	IO0	I/O	I/O
18	D3	VCC	VCC	DC-DC filtered output
19	D4	VBAT	VBAT	VBAT
20	D5	GND	GND	GND

## J102 SiP

- 4mm\*4mm\*0.75mm
- 20 pins
- MP: 2024. July




Pad #	Pad Name	Description	I/O	A/D
1	RF	RF signal input/output (Antenna)	I/O	
2	IO11	Digital input output 9	I/O	
3	IO9	Digital input output 10	I/O	
4	VCC	DC-DC filtered output		P
5	AGND	Analog ground		P
6	NRESET	Reset pin	I	
7	IO0	Digital input output / ADC 0	I/O	A/D
8	IO1	Digital input output / ADC 1	I/O	A/D
9	IO2	Digital input output / ADC 2	I/O	A/D
10	IO4	Digital input output 4	I/O	D
11	IO3	Digital input output / ADC 3	I/O	A/D
12	VBAT	Battery input voltage		P
13	IO10	Digital input output 10	I/O	D
14, 21	GND	Digital ground		P
15	IO12	Digital input output 12	I/O	D
16	IO13	Digital input output / CM3-JTAG Test Reset	I/O	D
17	TMS	CM3-JTAG Test Mode State	I/O	D
18	TCK	CM3-JTAG Test Clock	I/O	D
19	IO14	Digital input output / CM3-JTAG Test Data In	I/O	D
20	IO15	Digital input output / CM3-JTAG Test Data Out	I/O	D





# JHEAR J10/J102 Platform Introduction

Features	Comment
1-2-3 Buttons	Flexible choice for volume up/down
additional touch sensor interface	Customized touch function for volume, mode etc., better waterproof
50-0hm antenna port	Save PCB size and design cycle time
Support Multi customized Mode	In door, meeting, outdoor, Restaurant mode, music mode, ASHA mode
TV streamer Mode	Property protocol with ultra low power consumption, broadcasting, enable group flexibly
8-WDRC	8 channel WDRC
AI-base Noise reduction	AI NR model, will full operate at next generation platform
Feedback Algorithm	Support in pulse noise compression
ASHA function	Embedded ASHA function
Wireless fitting (App, Website)	reference code for App design, web remote fitting through <a href="https://yp.jhearing.com/webj10_2022.php">https://yp.jhearing.com/webj10_2022.php</a>
Rechargeable battery	Mature design reference, charging indicator
Pure Tone test	Open interface for self test
audio prompt function	Support custom multi-languages



# J10/J102 Wireless hearing aid Platform



## Basic Algorithm

- Customized 8-ch WDRC
- 16-Ch EQ
- Ultra low power audio streaming
- Starry AI Noise reduction
- AFC Feedback Control
- Howling Detect
- Pure-Tone Audiometry(PTA)

## Functions & Features

- TX SRC Sample Rate Convert: 16kHz
- 1->X style Audio broadcasting
- Marked Group for group education etc.
- Max latency < 6 ms
- TV streamer support 15 meters
- Rechargeable battery
- Remote Mic

## Solution support

- OTC hearing aid "turn-key" solution
- Open-source Mobile App
- PC Fitting source code sharing
- TV streamer(3.5mm, AUX, SPDIF)
- MFi support(special approval)
- ASHA function support
- Custom logo, name, audio prompt

# J10/J102 OTC Hearing Aid Scenarios

Support the coexistence of hearing aid features and ASHA function

- ❖ Volume control
- ❖ Mode switch
- ❖ EQ modification
- ❖ Pure Tone Test



Mobile Wireless Fitting



- ❖ WDRC
- ❖ NR
- ❖ Signal Chain Calibration
- ❖ Volume
- ❖ Mode
- ❖ AGC-O

PC Wireless Fitting

## J10 hearing aids



- ❖ TV-Box



- ❖ ASHA
- ❖ MFi



- ❖ TV-Box
- ❖ Remote Mic




# J10/J102 bases hearing aid use cases




## Direct connectivity with a Smartphone (no dongle)

- Streaming music and speech to earbuds/HA
- Controlling the earbuds/HA (volume up/down, equalizer etc..)



## Connectivity with accessories i.e. TV set (via a dongle)

- The earbuds/HA receive audio stream from TV
- The Smartphone controls the earbuds/HA




# J10/J102 Platform for Audio

## TV Steamer

With TV Steamer, you can

- watch TV, listen to music. we support audio broadcasting, one host for many slave device
- Support Group isolate group broadcasting
- Support all families TV enjoying together



## ASHA/MFi


Support most of Android phone, the system >10.0.0 version

- ASHA is audio streaming for Android smart phone
- This is for single direction audio receive




# J10 Platform Performance picture (OTC)


## ANSI specification to support OTC hearing aid benchmark




## Delay=6ms



## RIC hearing aid



## ITE hearing aid




# J10 OTC hearing aid design resource

<https://en.jhearing.com/docs/>

## J10 Document

- J10 Description
- Datasheet
- BLE protocol and client code
- J10 Schema pdf
- J10 Functions
- J10 Web fitting
- Downloads



[https://yp.jhearing.com/webj10\\_2022.php](https://yp.jhearing.com/webj10_2022.php)

**J10 Hearing Aid ,you can fitting on the web page**

Current Mode:

	Name	Band 0	Band 1	Band 2	Band 3	Band 4	Band 5	Band 6	Band 7
#	Cross over Frequency 频率		317Hz	502Hz	797Hz	1264Hz	2005Hz	3181Hz	5044Hz
#	Expansion Rate 放大系数	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>	<input type="text" value="1.0"/>
#	Expansion Knee point 放大区拐点	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>	<input type="text" value="30 dB"/>
#	Linear Gain 线性增益	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>	<input type="text" value="20 dB"/>
#	Compress Knee 压缩拐点	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>	<input type="text" value="55 dB"/>
#	Compress Ratio 压缩系数	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>	<input type="text" value="1.5"/>
#	Limit Knee 限制区拐点	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>	<input type="text" value="90 dB"/>

# J10 base power consumption

Hearing Aid work Status	Current	Comment		Summary
J10 Pure Hearing aid function	0.8mA	Include all algorithm functions	Feedback algorithm cost 0.7mA, better in next platform	<p>for this customer hearing aid product:</p> <ul style="list-style-type: none"> <li>❖ hearing aid algorithm total 0.8mA~,</li> <li>❖ DMIC cost 0.5mA~1mA</li> <li>❖ BLE cost 0.6~0.8mA</li> <li>❖ ASHA cost 2mA</li> </ul> <p>(Note: different customers do have different settings and output power for Audio streaming, the data could be for reference only)</p>
J10 support wireless ASHA	2.0mA	Feedback algorithm cost 0.7mA	We could reduce this when more powerful core in J20 platform by optimize the sample rate	
Shut down WDRC, EQ	2.5mA~ 2.6mA	The rest of hearing aid algorithm do not cost much power consumption		
System work, but shut down DMIC	2.3mA	DMIC power consumption in this hearing aid cost 1mA	There is better power consumption DMIC source in the future	
Only BLE, but no hearing aid algorithm	0.8mA	Which means if shut down BLE, we will save 0.8mA		



# J10 Hearing Aid Accessories

USB dongle



TV Streamer



TV Streamer



Remote Mic



3.5mm

SPDIF

HDMI IN

HDMI OUT

Support Ezairo 7160/J10 hearing aid

# Support JHEAR J10/ Ezairo 7160SL Hearing Aid Platform TV streamer







- High quality audio streaming from your TV, Hi-Fi stereo, personal computer or other audio device directly to your J10/Ezairo 7160 SL-based hearing aid
- TV streamer from JHEAR
  - Support HDMI, 3.5mm analog audio, coaxial, optical audio
  - Type-C charge interface, 5v-2A adapter
  - Support V+ / V- push button
  - Support Power on/off LED light and Audio signal indicator
- TV streamer can be customized if desired
  - Design, features and form factor can be modified
- To access the reference design:
  - JHEAR also offers full manufacturing services as well
  - JHEAR contact– [tech@jhearing.com](mailto:tech@jhearing.com)




# J10/J102 Hybrid SoC platform eco-system build for ramping up market

JHEAR could help support your HA supply chain to help your product best BOM and efficient design cycle



# JHEAR Hearing aid roadmap


- Released
- In development
- proposal




# J20 Over-The-Counter (OTC) Hearing Aids concept

## The Key content of J20

- Open Algorithm interface
  - Features customizable by customers
  - Kindly SDK
- Low-power Bluetooth Classic TWS earbud
  - A2DP with ear-to-ear relay
  - Hands-Free Profile support bi-directional audio phone calls
- LE Audio version for absolute lowest power
  - Isochronous stream to each ear (LC3 codec)
  - Auracast (broadcast) receiver
- Transparency mode
  - With EQ, noise reduction and directional processing
- Additional OTC hearing aid features
  - 16 channel WDRC
  - Feedback cancellation
  - AI base Noise reduction
  - Wireless fitting and remote control
  - Support natively Li-Ion and ZnAir battery
    - Li-Ion will use DCDC mode
    - ZnAir will use LDO mode



# J10 platform Success stories



More coming soon...




# JHEAR support whole supply Chain



## Antenna

### Antenna category

- LDS/PDS/LAP antennas
- Mobile phone/laptop/panel antenna
- TWS earphone/Hearing Aid antenna
- Internet of Things/Smart "Water and Electricity" Three Antennas
- NFC/RFID antennas
- V2X/UWB/GNSS antennas for new energy vehicles



## R&D Capability Resource

### 3D Chamber




Shenzhen 8 Set: ETS & SATIMO  
Chambers  
TW 1 set: ETS  
Shanghai 1set: ETS & Satimo  
XiAn 1 set: ETS  
Fuzhou 2 Set: ETS & Satimo

### Test equipment



Cellular protocol testing :BT, DVB-H, DVB-T and GPS, LTE SISO  
Positive test:  
2G,3G,4G,5G,WiFi/6E,BT,GPS  
NFC (Micropross / EMVCo)

### RF simulation Tool



CST MICROWAVE STUDIO  
HFSS 2 Sets  
ADS Circuit  
Maxwell

### Mechanical modeling software



PRO/E 5.0  
Catia  
Creo 8.0  
CAD 2019





Follow Us @ JHEAR



[www.jhearing.com](http://www.jhearing.com)

# J10 Algorithms flow

